

Expected medical instructions	The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger
Expected patient population	-People who need blood oxygen measurement .
Expected use or interaction with body parts tissue type	Finger
Expected user profile	People who need blood oxygen measurement , doctors, etc
Application environment	Avoid electromagnetic interference Extreme temperature Avoid pollution and dust Avoid direct sunlight, etc
Operating principle	Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of diodes through process in electronic circuits and microprocessor.

Classification

- 1.Management Class for Medical Devices: Class IIa
2.Anti-electric Shock Type: Internally powered equipment
3.Anti-electric Shock Degree: Type BF equipment
4.Overvoltage category classification:Class I
5.Pollution degree:
Pollution degree2:Micro-environment with non-conductive pollution,expect occasional conductivity caused by condensation

Maintenance and Preservation

- 1.Replace the batteries timely when the low battery indicator flashes
2.Clean the surface of fingertip pulse oximeter before use.
3.Remove the batteries when the oximeter is not likely to be used for some time when leakage from batteries would result in an unacceptable risk.
4.It would be better to preserve the product in 10-40°C (14-104°F) and humidity is 10%-95%.
5.It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
6.Cleaning frequency
If the oximeter is dirty when used at home, it is recommended to clean the enclosure and silicone pad after each use.If the oximeter is not dirty, simply clean the black silicone pad before and after each use.
When used in a medical institutions, clean it after each use.
7.Cleaning method
When the oximeter is used at home, wipe and disinfect with 75% alcohol, then dry naturally or clean the product with dry cloth. Clean at least twice a week.
When using in medical institutions, clean the tested finger with 75% alcohol before use. After each measurement, wipe and disinfect the contact part between black silicone and enclosure with 75% alcohol.
8.The oximeter can display functional arterial oxygen saturation and pulse rate after calibration.
9.The name of the simulator is FLUKE Index2 simulator, and the version number is 3.00.
10.The blood oxygen simulator is calibrated to display functional oxygen saturation.
11.The maximum temperature of the contact surface between the product and human body does not exceed 41 °C.

Waste dispoisa

- 1.Please follow local laws to dispose of waste scrap.
2. Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components, including used batteries and packaging box

Components

Components	Quantity
The host	1 set
Lanyard	1 pcs
User manual	1 pcs
AAA battery (optional)	2pcs

Manufacturer’ s Declaration of the EUT

- Statement:
1.The Fingertip Pulse Oximeter or user should use the product in the electromagnetic environment specified in the following table, otherwise it may cause abnormal operation of the product.
2. Fingertip Pulse Oximeter is a table-top equipment, it suitable for medical unit and home use.
3. Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm(12 inches)to any part of theoximeter, including cables specified by the manufacturer.Otherwise, degradation of the performance of this Fingertip Pulse Oximeter could result
4. Warning: Use of this Fingertip Pulse Oximeter adjacent to or stacked with other equipment should be avoided because it could result in improper observed to verify that they are operating normally.
5. Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Guidance and manufacturer’ s declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

1	Guidance and manufacturer’ s declaration – electromagnetic emission	
2	The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Fingertip Pulse Oximeter should assure that it is used in such an environment.	
3	Emissions test	Compliance
4	RF emissions CISPR 11	Group 1
5	RF emissions CISPR 11	Class B
6	Harmonic emissions IEC 61000-3-2	N/A
7	Voltage fluctuations /flicker emissions IEC 61000-3-3	N/A

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer’ s declaration – electromagnetic immunity		
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.		
Immunity test	EN 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15kV air
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (>95 % dip in UT) for 5 sec	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30A/m

Guidance and manufacturer’ s declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer’ s declaration – electromagnetic immunity		
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.		
Immunity	EN 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A

Table 9 - Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power(W)	Distance (m)	Immunity TEST LEVEL
385	380 -390	TETRA 400	Pulse modulation ^{b)}	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation	2	0.3	28
710	704 - 787	LTE Band 13,17	Pulse modulation ^{b)} 217 Hz	2	0.3	9
745						
780						
810	800 - 960	GSM 800/900,TETRA 800, iDEN 820,CDMA 850,LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1 170	1 700-1 990	GSM 1800;CDMA 1900; GSM 1900;DECT;LTE Band 1, 3 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400-2 570	Bluetooth,WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5 240	5 100-5 800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0.3	9
5 500						
5 785						

NOTE:

- If necessary to achieve the immunity test level, the distance between the transmitting antenna and the me equipment or me system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.
a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Possible Problems and Resolutions

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2.Patient’ s Oxyhemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unsteady	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status	1. Retry by plugging the finger 2. Please remain at rest
The Oximeter can not be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged	1. Please replace the batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected in 8 seconds 2. Inadequate power	1. Normal 2. Replace the batteries

Symbols and Definitions

	Type BF applied part		This device has no alarm function		separate collection for electrical and electronic equipment
	This way up		Humidity limitation		Not made with natural rubber latex
	Refer to instruction manual		Keep dry		First characteristic numeral 2: Against ingress of solid foreign objects: ≥12.5mm diameter Second characteristic numeral 2 Against ingree of water with harmful effects:dripping(15° tilted)
	Caution		Medical device		Model number
	Serial number		Batch code		Unique Device Identifier
	Date of manufacture		Temperature limit		signal inadequacy ① Indication of probe faults (open circuit condition or close circuit condition) ② Indication of Probe cable faults ③ Indication of Probe cable extender faults
	Manufacturer		Keep away from sunlight		
	Stand-by		Use-by date		

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

Statement:

1. If you need maintenance, please contact the manufacturer
2. The company can be in the form of email or other electronic files provide users with random files.
3. The instrument is not used for evaluation of blood oxygen probe pulse and pulse blood oxygen monitor accuracy.

After-sales service

Ensure that users

- Please read user manual before using the instrument;
- According to the requirement of the instruction manual for the operation and daily maintenance, and make sure the machine power supply, and environmental requirements

Maintenance regulations

- To conform to the regulations, free maintenance within the scope of products, with warranty card for free maintenance.All that is beyond the scope of free maintenance product, provide paid services.
- With warranty card and shopping invoice,main machine for a year, accessories for three months are under free maintenance services from the date of purchase.
- Following does not belong to the scope of free maintenance
 - ⊗ The fault caused by human factors, the damage;
 - ⊗ Due to the use to be inconsistent with the provisions of our company work environment to cause damage;
 - ⊗ Due to the product in the our company authorized personnel disassembling or repairing damaged;
 - ⊗ Products beyond the warranty period.